PRO powered by BRAUN

TRIZONE

OXYJET 1000-5000 SMARTSERIES

Type 3724  
Handle Type 3754

Deutsch	6
English	12
Français	17
Español	23
Português	29
Italiano	35
Nederlands	41
Dansk / Norsk	47
Svenska	53
Suomi	59
Ελληνικά	64

ⓃⓄ ⓃⓄ 00 800 27 28 64 63

ⓃⓂ 08 44 - 88 40 10

ⓄⓅ 0 800 731 1792

Ⓜⓔ 1 800 509 448

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Ⓞⓔ 0 800 14 592

ⓄⓁ 800 21172

ⓄⓈ 901 11 61 84

ⓄⓉ 808 20 00 33

ⓂⓉ 800 124 600

ⓄⓁ 0800-4 45 53 88

ⓄⓀ 70 15 00 13

ⓄⓃ 22 63 00 93

Ⓞⓔ 020 - 21 33 21

ⓄⓂ 020 377 877

ⓄⓇ 800 801 3457

**Internet:**

[www.oralb.com](http://www.oralb.com)  
[www.braun.com](http://www.braun.com)  
[www.service.braun.com](http://www.service.braun.com)

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# English

Welcome to Oral-B!

Please read the usage instructions carefully before first use and save this manual for future reference.

## IMPORTANT

- The appliance must be mounted in a way that it cannot fall into water or drop to the floor. The plug should always be accessible.
- Periodically check the cord for damage. If cord is damaged, take the charging unit to an Oral-B Service Centre. A damaged or non-functioning unit should no longer be used. Do not modify or repair the product. This may cause fire, electric shock or injury.
- The toothbrush is not intended for use by children under age three. Children aged from 3 years and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge can use the toothbrush, if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Cleaning and user maintenance shall not be performed by children.
- The oral irrigator is not intended for use by children under age 8. Children aged from 8 years and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge can use the product, if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children shall not play with the appliance.
- Oral irrigators shall be used only with liquids defined in the manufacturer's instructions.
- Use this product only for its intended use as described in this manual. Do

not use attachments or chargers which are not recommended by the manufacturer.

## WARNING

- If the product is dropped, the brush head/irrigator jet should be replaced before the next use even if no damage is visible.
- Do not place the appliance in water or liquid or store where it can fall or be pulled into a tub or sink. Do not reach for it when fallen into water. Unplug immediately.
- Do not disassemble the product except when disposing of the battery. When taking out the battery from the toothbrush for disposal of the unit, use caution not to short the positive (+) and negative (-) terminals.
- When unplugging, always hold the power plug instead of the cord. Do not touch the power plug with wet hands. This can cause an electric shock.
- If you are undergoing treatment for any oral care condition, consult your dental professional prior to use.
- This appliance is a personal care device and is not intended for use on multiple patients in a dental practice or institution.

## Important Information

- Your Oral-B SmartSeries toothbrush can be used with a Smart Phone (find details under «Connecting your toothbrush to your Smart Phone»).
- To avoid electromagnetic interference and/or compatibility conflicts, deactivate its radio transmission before using it in restricted surroundings, such as aircrafts or specially marked areas in hospitals. Deactivate radio transmission when your toothbrush is turned off by pressing the on/off button (5) for 3 seconds. All indicator lights (2, 3 & 6) will flash simultaneously. Follow the same procedure to reactivate transmission again. **Note:** Every time you turn your toothbrush off, when using it without radio transmission, the indicator lights will flash simultaneously.
- Persons with pacemakers should always keep the toothbrush more than six inches (15 centimeters) from the pacemaker while turned on. Anytime you suspect interference is taking place, deactivate the radio transmission of your toothbrush.

## Description

### Toothbrush

- 1 Charger base
- 2 Low Charge Indicator light
- 3 Charge Indicator light

- 4 Toothbrush handle (with *Bluetooth*® wireless technology)
- 5 On/off button (mode selection)
- 6 Pressure sensor light
- 7 Brush head
- 8 Brush head compartment

#### Irrigator

- 9 Jet compartment
- 10 Irrigator jet
- 11 Jet switch
- 12 Jet release button
- 13 On/off switch
- 14 Water pressure control dial
- 15 Handle switch (water stop and go)
- 16 Irrigator handle
- 17 Water container (600 ml)
- 18 Cord storage

#### Specifications

For voltage specifications please refer to the bottom of the appliance.

## A Using your toothbrush

### Connecting and charging

Your toothbrush has a waterproof handle, is electrically safe and designed for use in the bathroom.

- Plug the appliance into an electrical outlet and place the toothbrush handle (4) on the charge base (1) (a).
- The green charge indicator light (3) flashes while the handle is being charged. Once it is fully charged, the light turns off. A full charge can take up to 22 hours and enables up to seven days of regular brushing (twice a day, 2 minutes).  
Note: After deep discharge, the green light might not flash immediately; it can take up to 10–15 minutes.
- If the rechargeable battery is low the red low charge indicator light (2) is flashing for some seconds when turning your toothbrush on/off.
- For everyday use, the toothbrush handle can be stored on the Charger base to maintain it at full power. Overcharging is impossible.
- To maintain the maximum capacity of the rechargeable battery, unplug the charger base and fully discharge the handle with regular use at least every 6 months.

### Brushing techniques

Wet brush head and apply any kind of toothpaste. To avoid splashing, guide the brush head (7) to your teeth before switching on the appliance (b). When brushing your teeth with one of the **Oral-B oscillating-rotating brush heads** guide the brush head slowly from tooth to tooth, spending a few seconds on each tooth surface. When using the

**Oral-B «TriZone» brush head** place the toothbrush bristles against the teeth at a slight angle towards the gumline. Apply light pressure and start brushing in back and forth motions, just like you would do with a manual toothbrush. With any brush head start brushing the outsides, then the insides and finally the chewing surfaces. Brush all four quadrants of your mouth equally. Do not press too hard, simply let the brush do all the work. You may also consult your dentist or dental hygienist about the right technique for you.

In the first days of using any electric toothbrush, your gums may bleed slightly. In general, bleeding should stop after a few days. Should it persist after 2 weeks, please consult your dentist or dental hygienist. If you have sensitive teeth and gums, Oral-B recommends using the «Sensitive» mode (depending on model).

### Brushing modes (depending on model)

Model 1000	Model 2000	Model 3000	Model 4000/5000	Brushing modes (with handle type 3754)
✓	✓	✓	✓	<b>Daily Clean</b> – Standard for daily mouth cleaning
		✓	✓	<b>Sensitive</b> – Gentle, yet thorough cleaning for sensitive areas
			✓	<b>Whitening</b> – Exceptional polishing for occasional or everyday use
	✓	✓	✓	<b>Gum Care</b> – Gentle massage of gums

Your toothbrush automatically starts in the «Daily Clean» mode. To switch to other modes, successively press the on/off button (5). If you wish to turn off your toothbrush from «Sensitive» mode e.g., push and hold down the on/off button until the motor stops.

### Professional Timer

A short stuttering sound at 30-second intervals reminds you to brush all four quadrants of your mouth equally (d). A long stuttering sound indicates the end of the professionally recommended 2-minute brushing time. The elapsed brushing time is memorised even when the handle is briefly switched off during brushing. When pausing longer than 30 seconds the timer resets.

### Pressure sensor

Protecting your teeth and gums from hard brushing, your toothbrush has a pressure control feature installed. If too much pressure is applied, the red pressure sensor light (6) will light up reminding you to reduce pressure (c). In addition you will also hear

a different sound while brushing. Periodically check the operation of the pressure sensor by pressing lightly on the brush head during use.

### Brush heads

Oral-B offers you a variety of different brush heads that fit your Oral-B toothbrush handle.

Our **oscillating-rotating brush heads** can be used for tooth-by-tooth precision cleaning.

#### Oral-B CrossAction brush head

Our most advanced brush head. Angled bristles for a precise clean. Lifts and powers away plaque.

#### Oral-B FlossAction brush head

features micropulse bristles that allow for an outstanding plaque removal in the interdental areas.

#### Oral-B 3D White brush head

features a special polishing cup to naturally whiten your teeth. Please note that children under 12 years old should not use the Oral-B «3D White» brush head.

#### Oral-B Sensitive brush head

features an extra soft bristle structure that is gentle on teeth and gums.

#### Oral-B PowerTip brush head

offers a deep clean in hard to reach areas, e.g. around bridges, crowns or implants.

Our **sweeping «Trizone» brush head** can be used with all common brushing techniques.

#### Oral-B TriZone brush head

features triple zone cleaning action for outstanding plaque removal, even in between teeth.

Oral-B brush heads feature light blue INDICATOR® bristles to help you monitor brush head replacement need. With thorough brushing, twice a day for two minutes, the blue color will fade halfway approximately within 3 months, indicating the need to replace your brush head. If the bristles splay before the color recedes, you may exert too much pressure on teeth and gums.

We do not recommend using the Oral-B FlossAction or Oral-B 3D White brush head with braces. You may

use the Oral-B Ortho brush head, specifically designed to clean around braces and wires.

### Connecting your toothbrush to your Smart Phone

Download the Oral-B™ app for free from App Store<sup>(SM)</sup> or Google Play™. The Oral-B™ app gives you the opportunity to use your time more efficiently, chart your brushing progress and customize your toothbrush settings (for more benefits see the Oral-B™ app instructions).

- Before using the Oral-B™ app make sure *Bluetooth* is activated on your smart phone (for directions please refer to your smart phone's user manual). Your toothbrush will appear in the list of devices as: «Oral-B toothbrush». Note: Oral-B™ app functionality is limited if *Bluetooth* is deactivated.
- Start the Oral-B™ app.
- Once you start using your toothbrush by lifting it off the plugged-in charger base (1) or briefly turning your toothbrush on and off (5), it automatically connects to the Oral-B™ app.
- Any further app instructions will be displayed on your smart phone.
- Keep your smart phone in close proximity when using it with your toothbrush handle. Make sure the smart phone is safely positioned in a dry place.  
Note: Your smart phone must support *Bluetooth* 4.0 (or higher) in order to connect to your toothbrush handle.

### Cleaning recommendations

After brushing, rinse your brush head carefully under running water for some seconds with the handle switched on. Switch off the handle and remove the brush head. Clean both parts separately under running water (e); then wipe them dry before reassembling the toothbrush. Put the toothbrush handle back on the charger base. The charger base should be cleaned with a damp cloth only.

## B Using your irrigator

Before using the irrigator first time, fill the water container (17) and allow the water to run through completely by holding the irrigator handle (16) over a washbasin.

- Attach one of the irrigator jets (10) to the handle until it locks into place (f). Fill the water container with lukewarm water. If you want to add mouthwash into the water, please make sure you fill the water container with lukewarm water first, then add the mouthwash (g).
- Turn on the irrigator by pressing the on/off switch (13) (h). The handle switch should be in the lower position = water stop.

- Select the water pressure with the pressure control dial (14) (1 = gentle, 5 = strong) (i). We recommend that you start with setting 1, especially if your gums tend to bleed or if children are using the irrigator.
- Bend over a washbasin and guide the jet to your teeth and gums. Slide the handle switch (15) up to start the water stream (j). Keep your mouth slightly open to allow the water to run into the washbasin.

You can choose between two different jet functions by sliding the jet switch (11) down or up (before doing so, turn off the water stream by sliding the handle switch down):

Rotating stream for general cleaning and gum massage (jet switch in lower position): The rotation of the water stream and its micro-bubbles is designed for easy penetration and gentle gum massage (l).

Mono stream (jet switch in upper position): The straight water stream is designed for targeted cleaning (m).

Maximal operation time: 15 minutes

Cooling off time: 2 hours

In general we recommend to replace the irrigator jet every 3 months.

### After use

Always empty the water container entirely. If your dentist has recommended the use of mouthwash with your irrigator, allow some clean water to run through the irrigator before you empty the container entirely to prevent any deposit build-up. Make sure you keep the irrigator switched on until the water stops running through the irrigator jet, to prevent bacteria from collecting in the standing water.

Press the on/off switch to turn off the irrigator.

To take off the irrigator jet, press the jet release button (12) (k). Always wipe the handle and the irrigator jet dry; place the handle back onto its base and store the irrigator jet in the jet compartment (9).

### Cleaning recommendations

Regularly clean the housing with a damp cloth. The water container and the jet compartment are removable and dishwasher-safe.

Subject to change without notice.

## Environmental Notice

Product contains batteries and/or recyclable electric waste. For environment protection do not dispose in household waste, but for recycling take to electric waste collection points provided in your country.

## Guarantee

We grant 2 years guarantee on the product commencing on the date of purchase. Within the guarantee period we will eliminate, free of charge, any defects in the appliance resulting from faults in materials or workmanship, either by repairing or replacing the complete appliance as we may choose.

This guarantee extends to every country where this appliance is supplied by Braun/Oral-B or its appointed distributor.

This guarantee does not cover: damage due to improper use, normal wear or use, especially regarding the brush heads, as well as defects that have a negligible effect on the value or operation of the appliance. The guarantee becomes void if repairs are undertaken by unauthorised persons and if original Braun parts are not used.

This appliance is equipped with a *Bluetooth* Smart approved radio module. Failure to establish a *Bluetooth* Smart connection to specific smart phones is not covered under the device guarantee, unless the *Bluetooth* radio module of this appliance is damaged.

The *Bluetooth* devices are warranted by their own manufacturers and not by Oral-B. Oral-B does not influence or make any recommendations to device manufacturers, and therefore Oral-B does not assume any responsibility for the number of compatible devices with our *Bluetooth* systems. Oral-B reserves the right, without any prior notice, to make any technical modifications or changes to system implementation of device features, interface, and menu structure changes that are deemed necessary to ensure that Oral-B systems function reliably.

To obtain service within the guarantee period, hand in or send the complete appliance with your sales receipt to an authorised Oral-B Braun Customer Service Centre.

### For UK only:

This guarantee in no way affects your rights under statutory law.

### Important information on *Bluetooth*<sup>®</sup> radio module

Although all indicated functions on the *Bluetooth* device are supported, Oral-B does not ensure 100% reliability in connection, and feature operation consistency. Operation performance and connection reliability are direct consequences of each individual *Bluetooth* device, software version, as well as the operating system of those *Bluetooth* devices and company security regulations implemented on the device.

Oral-B adheres to and strictly implements the *Bluetooth* standard by which *Bluetooth* devices can communicate and function in Oral-B toothbrushes. However, if device manufacturers fail to implement this standard, *Bluetooth* compatibility and features will be further compromised and the user may experience functional and feature-related issues. Please note that the software on the *Bluetooth* device can significantly affect compatibility and operation.

**Warranty replacement brush heads**

Oral-B warranty will be voided if electric rechargeable handle damage is found to be attributed to the use of non-Oral-B replacement brush heads.

Oral-B does not recommend the use of non Oral-B replacement brush heads.

- Oral-B has no control over the quality of non-Oral-B replacement brush heads. Therefore, we cannot ensure the cleaning performance of non-Oral-B replacement brush heads, as

communicated with the electric rechargeable handle at time of initial purchase.

- Oral-B cannot ensure a good fit of non-Oral-B replacement brush heads.
- Oral-B cannot predict the long-term effect of non-Oral-B replacement brush heads on handle wear.

All Oral-B replacement brush heads carry the Oral-B logo and meet Oral-B’s high quality standards. Oral-B does not sell replacement brush heads or handle parts under any other brand name.

**30-Days Money Back Guarantee**

Take the Oral-B 30 day Challenge: Satisfaction guaranteed or your money back!

Full details at [www.oralb.com](http://www.oralb.com) or call:

**GB** 0800 731 1792

**IE** 1800 509 448

**Trouble Shooting**

Problem	Possible reason	Remedy
<b>TOOTHBRUSH</b>		
Gums bleed slightly	<ol style="list-style-type: none"> <li>1. May happen in first days of usage</li> <li>2. Brush head</li> <li>3. Brushing mode (depending on model)</li> </ol>	<ol style="list-style-type: none"> <li>1. Bleeding should stop after a few days. Should it persist after 2 weeks, please consult your dentist or dental hygienist.</li> <li>2. Try using the Oral-B «Sensitive» brush head.</li> <li>3. Try using the «Sensitive» mode (depending on model).</li> </ol>
<b>BRUSH HEADS</b>		
Oral-B brush head does not fit Oral-B handle.	Oral-B brush head not compatible with your handle.	Oral-B brush heads that fit your Oral-B handle do not fit «Pulsonic», «VitalitySonic», «Sonic Complete», «CrossAction Power».
<b>APP</b>		
Toothbrush does not operate (properly) with Oral-B app.	<ol style="list-style-type: none"> <li>1. Oral-B app is turned off.</li> <li>2. <i>Bluetooth</i> is not activated on your smart phone.</li> <li>3. Radio transmission is deactivated.</li> <li>4. Connection to <i>Bluetooth</i> has been lost.</li> <li>5. Your smart phone does not support <i>Bluetooth</i> 4.0 (or higher)</li> </ol>	<ol style="list-style-type: none"> <li>1. Turn on Oral-B app.</li> <li>2. Activate <i>Bluetooth</i> on your smart phone; described in its user manual.</li> <li>3. Activate the radio transmission by pressing the on/off (5) for 3 seconds.</li> <li>4. Briefly turn your toothbrush on and off (5) to re-connect. Keep your smart phone in close proximity when using it with the toothbrush handle.</li> <li>5. Your smart phone must support <i>Bluetooth</i> 4.0 (or higher) in order to connect to your toothbrush handle.</li> </ol>
Reset to factory settings.	Original functionality wanted.	Press and hold on/off button (5) for 10 seconds until all indicator lights (2, 3 & 6) will double blink simultaneously. Note: radio transmission may be (de)activated while resetting. After reset radio transmission is re-activated.