## English

- 1 Insert batteries as shown. Uses two 1.5 V type AA (LR 6). For best performance, we recommend Duracell® Ultra alkaline batteries. Batteries will last approx. 3 months when brushing twice a day for two minutes.
- 2 Remove protection cap, wet brushhead and apply toothpaste.
- 3 Place brushhead on teeth and switch on.
- 4 Brush for at least 2 minutes (recommended by dentists). Guide the brushhead from tooth to tooth, first the outsides (a), then the insides (b), finally the chewing surfaces (c).  
During the first days of use, gums may bleed slightly. Should bleeding persist after 2 weeks, contact your dentist.
- 5 Replace brushhead when indicated (a = new, b/c = replace)
- 6 Clean as shown under running water, then dry.
- 7 Dispose of discharged batteries according to national or local regulations.

Adult supervision recommended when used by children.

Subject to change without notice.

This product conforms to the EMC-Requirements as laid down by the Council Directive 89/336/EEC.

Manufacturer's guarantee: In case of faults in materials or workmanship, the appliance will be replaced within 2 years as of date of purchase. Not covered: use and wear of batteries and brushhead.

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