Oral-B®

Vitality™

D 12523

D 12513

D 12013

powered by **BRAUN**

Type 3709

(D) (A) 00 800 27 28 64 63

(CH) 08 44 - 88 40 10

(GB) 0 800 73 11 792

(IRL) 1 800 509 448

(F) 0 810 365 855

(B) 0 800 14 592

(LU) 4 00 50 51

(E) 901 11 61 84

(P) 808 20 00 33

(I) 02-667 86 23

(NL) 0800-4 45 53 88

(GR) 210 9478700

Internet:
www.oralb.com
www.braun.com

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English

Please read the use instructions carefully before first use.

Important

- Periodically check the cord for damage. A damaged or non-functioning unit shall no longer be used. Take the charging unit (E) to an Oral-B Braun Service Centre.
- This product is not intended for use by children under age three.
- If the product is dropped, the brushhead should be replaced before the next use even if no damage is visible.

Description

- A Brushhead (depending on model) – interchangeable
- B On switch
- C Off switch
- D Handle
- E Charging unit
- F Brushhead holder (D 12523 model only)

Charging

- Plug the charging unit (E) into an electrical outlet. Place the handle on the charging unit. Handle must be switched off.
- A full charge takes 16 hours and will provide an operation time of approximately 20 minutes.

- For everyday use, the handle can be stored on the plugged-in charging unit to maintain the toothbrush at full power. Overcharging is impossible.

Battery maintenance

To maintain the maximum capacity of the rechargeable battery, you must discharge it at least once every 6 months. To do this, unplug the charging unit and allow the battery to run down with normal use.

Using the toothbrush

- Wet brushhead and apply toothpaste.
- Place brushhead on teeth, then switch on.
- Slowly guide the brushhead from tooth to tooth, first the outsides, then the insides, finally the chewing surfaces. Do not press too hard or scrub.
- After two minutes (brushing time recommended by dentists) the timer (D 12523, D 12513 model only) will signal by a short stutter.
- During the first days of use, gums may bleed slightly. Should bleeding persist after 2 weeks, contact your dentist.
- Replace brushhead every 3 months, when bristles are worn (1c) or when colour on Indicator® bristles (1b) has faded half-way down the bristles.

Cleaning

After use, rinse the brushhead under running water (2). Remove the brushhead and rinse brushhead and handle separately (3), then dry.

Subject to change without notice.

This product conforms to the European Directives EMC 89/336/EEC and Low Voltage 2006/95/EC.

This product contains rechargeable batteries. In the interest of protecting the environment, please do not dispose of the product in the household waste at the end of its useful life. Disposal can take place at an Oral-B Braun Service Centre or at appropriate collection points provided in your country.

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